

## JOEY CLAMP™

The Joey Clamp™ incorporates several technologies that are hidden in its jewelry-like appearance. The sinusoidal waveform of the closure allows the clamp to fit automatically on virtually any umbilical cord diameter and to other physiological changes that occur following cord transection. The hinge and clasp are cleverly incorporated in the marsupial's ears. The Joey Clamp™'s gentle, rounded design reassures the mother with its friendly shape.



- Attractive Design
- Streamlined Clasp
- Integral Hinge
- Latched Closure
- Cute Baby Marsupial Face
- Colorful Polycarbonate Material
- Friendly Round Shape

## JOEY CLAMP & CUTTER™

The Joey Clamp & Cutter™ includes the Joey Clamp™, an umbilical cord cutter and a maternal clamp. The Joey Clamp & Cutter™ is ergonomically designed for one-handed, ambidextrous use. Its double latching system ensures a complete cut and closure of the device. The cutter barrel reduces blood spray, increasing the safety of the delivery team. Once fully latched, simply pressing outward on the marsupial's ears easily ejects the Joey Clamp™.



- Joey Clamp™, Cutter and Maternal Clamp
- One-handed, Ambidextrous Use
- Textured Surface
- Enclosed Cutting Space
- Translucent Polycarbonate Material
- Double-latching System
- Reduces Blood Spray

## ORDER INFORMATION

Part Number	Description	Color	Size	Quantity
JC-10-01	Joey Clamp™ - Sterile	Purple	Universal	125 ct./case
JCC-20-01	Joey Clamp & Cutter™ - Sterile	Purple	Universal	60 ct./case



The Utopia Group LLC is the exclusive licensee of the Joey Clamp and Cutter  
- [www.joeymedical.com](http://www.joeymedical.com)  
513.791.1910